


Vital Signs Monitor Operation Manual



— English
Version 1.1, 2017
All rights reserved.

Table of Contents

Table of Contents.....	I
Warranty and Service Information.....	V
Proprietary Notice.....	V
Limited Warranty.....	V
Service Support.....	V
Chapter 1: Introduction.....	1-1
1.1 About this Manual.....	1-1
1.2 Definition of Symbols.....	1-1
1.3 Warning Information.....	1-2
Chapter 2: Intended Use and General Information.....	2-1
2.1 Intended Use.....	2-1
2.2 Monitor Features.....	2-1
2.3 Theory of Operation.....	2-2
2.3.1 Oxygen Saturation (SpO2).....	2-2
2.3.2 End-Tidal Carbon Dioxide Concentration (ETCO2).....	2-3
Chapter 3: Controls and Features.....	3-1
3.1 Monitor Front Panel.....	3-1
3.2 Monitor Operating keys.....	3-3
Chapter 4: Operating Instructions.....	4-1
4.1 Unpacking the Monitor.....	4-1
4.2 Install the Batteries.....	4-1
4.3 AC Power Adapter.....	4-2
4.3.1 Charging Ni-MH batteries.....	4-2
4.4 Attaching the Sensor to the Patient.....	4-2
4.4.1 Choosing the Sensor.....	4-3
4.4.2 Care and Handling of the Sensor.....	4-3
4.4.3 Checking the Sensor and Adapter Cable.....	4-7
4.4.4 Cleaning or Disinfecting the Sensors.....	4-9
4.5 Performance Considerations.....	4-9
4.5.1 SpO2 Performance Considerations.....	4-9
4.5.2 ETCO2 Performance Considerations.....	4-10
4.6 Turning On the Monitor.....	4-10
4.6.1 Six Display Modes.....	4-12

4.6.2 Two Measuring Modes.....	4-13
4.7 Turning Off the Monitor.....	4-13
Chapter 5: Changing the Monitor's Settings.....	5-1
5.1 Changing the Pulse Volume.....	5-1
5.2 Changing the Alarm Volume.....	5-1
5.3 Managing the Patient's Information.....	5-1
5.4 Changing the Alarm Limits.....	5-2
5.4.1 Alarm Summary.....	5-2
5.4.2 Changing the SpO2 Alarm Settings.....	5-3
5.4.3 Changing the ETCO2 Alarm Settings.....	5-4
5.5 Changing ETCO2 Settings.....	5-5
5.6 System Settings.....	5-7
5.6.1 Changing the ETCO2 Alarm Settings.....	5-7
5.7 Resetting.....	5-8
Chapter 6: Trend Data Output and Analysis.....	6-1
6.1 Connecting the Device to PC.....	6-1
6.2 Installing the Data Management System.....	6-1
6.3 Function of Main Menu.....	6-2
6.4 Trend Data Download.....	6-2
6.5 Open the Data Base.....	6-4
6.6 Data Analysis and Reports.....	6-5
Chapter 7: Accessories.....	7-1
7.1 Standard Configuration.....	7-1
7.2 Optional Accessories.....	7-1
Chapter 8: Maintenance and Troubleshooting.....	8-1
8.1 Schedule of Maintenance.....	8-1
8.2 Storage.....	8-1
8.3 Troubleshooting.....	8-2
8.4 EMC (Electro-Magnetic Compatibility).....	8-4
Chapter 9: Specification.....	9-1
9.1 Equipment Classification.....	9-1
9.2 Parameter.....	9-1
9.6 Default Settings of Alarms Limits.....	9-1
9.7 Power Requirements.....	9-1

9.8 Battery Life.....	9-2
9.9 Dimensions.....	9-2
9.10 Environmental Specification.....	9-2

Warranty and Service Information

Proprietary Notice

Information contained in this document is copyrighted by the manufacturer, Ltd. and may not be duplicated in full or part by any person without prior written approval of manufacturer. Its purpose is to provide the user with adequately detailed documentation to efficiently install, operate, maintain and order spare parts for the device supplied. All information contained in this document is believed to be current and accurate as of the date of publication or revision, but does not constitute a warranty.

Limited Warranty

The manufacturer ("Seller") warrants each new device to be free from defects in workmanship and materials under normal use and service for a period of one (1) years from the date of shipment. The sole obligation of manufacturer under this warranty will be repair or replace, at its option, products that prove to be defective during the warranty period. The foregoing shall be the sole warranty remedy. Except as set forth herein, seller makes no warranties, either expressed or implied, including the implied warranties of merchantability and fitness for a particular purpose. No warranty is provide if the products are modified without the express written consent of manufacturer and seller shall not be liable in any event for incidental or consequential damage. This warranty is not assignable.

Warranties are subject to change. Please contact manufacturer for current warranty information.

Service Support

Repairs for devices manufactured by manufacturer under warranty must be made at authorized repair centers. If the device needs repair, contact your local distributor or the manufacturer after-service department. When calling, have the device's model and serial number ready.

If you need to ship the device, pack the device and accessories carefully to prevent shipping damage. All accessories should accompany the device.

NOTE! Shipments received without a return number will be returned to sender.













Chapter 1: Introduction

1.1 About this Manual

The Operation Manual provides installation, operation, and maintenance instructions for health-care professionals and other users, trained in monitoring respiratory and cardiovascular activity.

These instructions contain important information for the safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Cautions, before using the monitor. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient.

1.2 Definition of Symbols

SYMBOLS	DEFINITION
	Attention, see in instructions for use
	Type BF Defibrillation
	Power on/off
	Alarm silence
	Up and Down Arrows
	Mode Change Key
	Menu Key
	Date of Manufacturing
IPX1	Drip Proof (monitor only)
	Indicates separate collection for electrical and electronic equipment.
	Probe is not connected
	Finger and sensor are not properly aligned
	Searching pulse
CO2 X	CO2 module has not been activated.
CO2 ✓	CO2 module has been activated.

1.3 Warning Information

KEYWORD	DEFINITION
WARNING	WARNING indicates a potential hazard or unsafe practice which, if not avoided, could result in death or serious injury.
CAUTION	CAUTION indicates a potential hazard or unsafe practice which, if not avoided, could result in minor personal injury or product /property damage.
NOTE	NOTE provides application tips or other useful information to assure that you get the most from your equipment.

Warnings

WARNING! Do not use this device in the presence of flammable anesthetics or other flammable substance in combination with air, oxygen-enriched environments, or nitrous oxide.

WARNING! Do not use this device in the presence of magnetic resonance imaging (MR or MRI) equipment.

WARNING! Use only SpO₂ sensors supplied with, or specifically intended for use with, this device.

WARNING! Do not plug the monitor into an outlet controlled by a wall switch.

WARNING! This device is intended for use by persons trained in professional health care. The operator must be thoroughly familiar with the information in this manual before using the device.

WARNING! This device must be used in conjunction with clinical signs and symptoms. This device is only intended to be an adjunct in patient assessment.

WARNING! It is the operator's responsibility to set alarm limits appropriately for each individual patient.

WARNING! Prolonged use or the patient's condition may require changing the sensor site periodically. Change sensor site and check skin integrity, circulatory status, and correct alignment at least every 4 hours.

WARNING! ELECTRICAL SHOCK HAZARD when cover is removed. Do not

remove covers. Refer servicing to qualified personnel.

WARNING! Chemicals from a broken LCD display panel are toxic when ingested. Use caution when handling a monitor with a broken display panel.

WARNING! In the event that earth ground integrity is lost, the performance of this device and/or other devices nearby may be affected due to excessive RF emissions.

WARNING! Any monitor that has been dropped or damaged should be inspected by qualified service personnel, prior to use, to insure proper operation.

WARNING! If the accuracy of any measurement is in question, verify the patient's vital sign(s) by an alternative method and then check the monitor for proper functioning.

WARNING! Remove device batteries prior to long term storage.

WARNING! Operation of this device may be adversely affected in the presence of strong electromagnetic sources, such as electrosurgery equipment.

WARNING! Operation of this device may be adversely affected in the presence of computed tomograph (CT) equipment.

WARNING! Use only SpO₂ sensors supplied with, or specifically intended for use with, this device.

WARNING! SpO₂ measurements may be adversely affected in the presence of high ambient light. Shield the sensor area (with a surgical towel, for example) if necessary.

WARNING! Dyes introduced into the bloodstream, such as methylene blue, indocyanine green, indigo carmine, patent blue V (PBV), and fluorescein may adversely affect the accuracy of the SpO₂ reading.

WARNING! Any condition that restricts blood flow, such as use of a blood pressure cuff or extremes in systemic vascular resistance, may cause an inability to determine accurate pulse rate and SpO₂ readings.

- WARNING!** Significant levels of dysfunctional hemoglobins, such as carboxyhemoglobin or methemoglobin, will affect the accuracy of the SpO₂ measurement.
- WARNING!** The monitor was not designed or tested to be an apnea monitor.
- WARNING!** Optical cross-talk can occur when two or more sensors are placed in close proximity. It can be eliminated by covering each site with an opaque material.
- WARNING!** Tissue damage may result from overexposure to sensor light during photodynamic therapy with agents such as verteporphin, porfimer sodium, and metatetrahydroxyphenylchlorin (mTHPC). Change the sensor site at least every hour and observe for signs of tissue damage. More frequent sensor site changes/inspections may be indicated depending upon the photodynamic agent used, agent dose, skin condition, total exposure time or other factors. Use multiple sensor sites.
- WARNING!** When connecting this monitor to any instrument, verify proper operation before clinical use. Refer to the instrument's user manual for full instructions. Accessory equipment connected to the monitor's data interface must be certified according to the respective IEC standards, i.e., IEC 60950 for data processing equipment or IEC 601-1 for electromedical equipment. All combinations of equipment must be in compliance with IEC 601-1-1 systems requirements. Anyone connecting additional equipment to the signal input port or the signal output port configures a medical system, and, therefore, is responsible that the system complies with the requirements of the system standard IEC 601-1-1.
- WARNING!** Disconnect the AC power supply from the outlet before disconnecting it from the monitor. Leaving the AC power supply connected to an AC power outlet without being connected to the monitor may result in a safety hazard.
- WARNING!** Do not allow any moisture to touch the AC power supply connectors or a safety hazard may result. Ensure that hands are thoroughly dry before handling the AC power supply.
- WARNING!** Failure to place the monitor away from the patient may allow the

patient to turn off, reset, or damage the monitor, possibly resulting in the patient not being monitored.

- WARNING!** Failure to carefully route the cable from the sensor to the monitor may allow the patient to become entangled in the cable, possibly resulting in patient strangulation. Route the cable in a way that will prevent the patient from becoming entangled in the cable. If necessary, use tape to secure the cable.
- WARNING!** If there is a risk of the AC power supply becoming disconnected from the monitor during use, secure the cord to the monitor several inches from the connection.
- WARNING!** Under certain clinical conditions, the device may display dashes if unable to display SpO₂ and/or pulse rate and/or ETCO₂ values. Under these conditions, the device may also display erroneous values. These conditions include, but are not limited to: patient motion, low perfusion, cardiac arrhythmias, high or low pulse rates or a combination of the above conditions. Failure of the clinician to recognize the effects of these conditions on the device readings may result in patient injury.
- WARNING!** Verify that all LEDs (light emitting diodes) on the display light up upon startup of the device.
- WARNING!** Explosion Hazard: DO NOT use in the presence of flammable anesthetics. Use of the CO₂ Sensor in such an environment may present an explosion hazard.
- WARNING!** Electrical Shock Hazard: Always disconnect the CO₂ Sensor before cleaning. DO NOT use if it appears to have been damaged. Refer servicing to qualified service personnel.
- WARNING!** Electrical Shock Hazard: Always disconnect the CO₂ Sensor Follow precautions for electrostatic discharge (ESD) and electromagnetic interference(EMI) to and from other equipment.
- WARNING!** Failure of Operation: If the CO₂ Sensor fails to respond as described in this userguide; DO NOT use it until approved for use by qualified personnel.
- WARNING!** Reuse, disassembly, cleaning, disinfecting or sterilizing the single

patient use CO2 airway adapters may compromise functionality and system performance leading to a user or patient hazard. Performance is not guaranteed if an item labeled as single patient use is reused.

WARNING! Inspect the CO2 airway adapters for damage prior to use. **DO NOT** use the CO2 airway adapters if they appear to be damaged or broken.

WARNING! Replace the CO2 airway adapters if excessive secretions are observed.

WARNING! If the CO2 waveform (Capnogram) appears abnormal, inspect the CO2 airway adapters and replace if needed.

WARNING! Monitor the CO2 waveform (Capnogram) for elevated baseline. Elevated baseline can be caused by sensor or patient problems.

WARNING! Periodically check the CO2 sensor and tubing for excessive moisture or secretion buildup.

WARNING! While using the CO2 sensor, a system leak, such as that caused by an uncuffed endotracheal tube or a damaged CO2 sensor may significantly affect flow-related readings. These include flow, volume, pressure, and other respiratory parameters.

WARNING! Do not operate the CO2 Sensor when it is wet or has exterior condensation.

WARNING! Inspect the gas dryer line and nasal sampling line for damage prior to use. **DO NOT** use the CO2 gas dryer line and nasal sampling line if they appear to be damaged or broken.

WARNING! Periodically check the CO2 sensor, gas dryer line and nasal sampling line for excessive moisture or secretion buildup.

Cautions

CAUTION! Do not autoclave, ethylene oxide sterilize, or immerse the monitor or sensors in liquid. Always disconnect the power source and remove all batteries before cleaning or disinfecting the monitor.

- CAUTION!** Should the device become wet, wipe off all moisture and allow sufficient time for drying before operating.
- CAUTION!** The monitor should be operated from its internal power source if the integrity of the protective earth conductor is in doubt.
- CAUTION!** Pressing front panel keys with sharp or pointed instruments may permanently damage the keypad. Press front panel keys only with your finger.
- CAUTION!** Do not allow water or any other liquid to spill onto the monitor. Unplug the external power supply from the monitor before cleaning or disinfecting the monitor. Evidence that liquid has been allowed to enter the monitor voids the warranty.
- CAUTION!** Ensure the device's AC rating is correct for the AC voltage at your installation site before using the monitor. The monitor's AC rating is shown on the external power supply. If the rating is not correct, do not use the monitor; contact the manufacturer after-service department, or your local distributor, for help.
- CAUTION!** Chemicals used in some cleaning agents may cause brittleness of plastic parts. Follow cleaning instructions in this manual.
- CAUTION!** Electrical Shock Hazard; the CO₂ Sensor is not user serviceable.
- CAUTION!** Follow precautions for electrostatic discharge (ESD) and electromagnetic interference (EMI) to and from other equipment.
- CAUTION!** Law restricts the sale, distribution, or use of CO₂ Sensors to or by the order of a licensed medical practitioner.
- CAUTION!** Use only National Medical approved accessories
- CAUTION!** DO NOT use the CO₂ Sensor if it appears to be damaged. Refer servicing to qualified personnel.
- CAUTION!** DO NOT use the CO₂ Sensor if it fails to properly operate.
- CAUTION!** DO NOT sterilize or immerse the CO₂ Sensor in liquids.

- CAUTION!** DO NOT clean the CO2 Sensor except as directed in this guide.
- CAUTION!** Avoid possible damage to the CO2 Sensor by following the cleaning and disinfection instructions in this guide.
- CAUTION!** Do not apply excessive tension to the CO2 Sensor cable.
- CAUTION!** Do not apply excessive tension to any sensor cable or pneumatic tubing.
- CAUTION!** Excessive moisture in the CO2 sensor may affect the accuracy of the flow measurement.
- CAUTION!** To avoid the affects of excessive moisture in the measurement circuit, insert the CO2 sensor in the ventiator circuit with the tubes upright.
- CAUTION!** It is recommended that the CO2 sensor be removed from the circuit whenever an aerosolized medication is delivered. This is due to the increased viscosity of the medications which may contaminate the sensor windows, causing the sensor to fail prematurely.
- CAUTION!** Do not store the CO2 Sensor at temperatures less than -40°F (-40°C) or greater than 131°F(55°C).
- CAUTION!** Do not operate the CO2 Sensor at temperatures less than 32°F (0°C) or greater than 113°F(45°C).
- CAUTION!** The Gas Dryer line can be used only 150 hours. If use it too long, it may cause wrong data or damage the module.
- CAUTION!** Do not use the module without connecting with the Drying line. Otherwise, it will cause measurement error or damage the module.

Notes

- NOTE!** Batteries are user replaceable. Follow local governing ordinances and recycling instructions regarding disposal or recycling of device components, including batteries.
- NOTE!** When using AC power, the device is a class II device with functional earth. This earth connection is for device electromagnetic compatibility

and does not provide protection to the patient or user.

NOTE! It is recommended that batteries be used with the monitor when operating with AC power to prevent monitor shutdown with loss of AC power.

NOTE! All user and patient accessible materials are non-toxic.

NOTE! Each input and output connection of the monitor is electrically isolated.

NOTE! The CO₂ Sensor contains no user serviceable parts. Refer service to qualified service personnel.

NOTE! Components of this product and its associated accessories which have patient contact are free of latex.

NOTE! Disposal of the CO₂ Sensor and its accessories should comply with national and/or local requirements.

NOTE! Nitrous oxide, elevated levels of oxygen, helium and halogenated hydrocarbons can influence the CO₂ measurement.

NOTE! Barometric pressure compensation is required to meet the stated accuracy of the CO₂ Sensor.

NOTE! As with all flow measuring devices, adverse conditions may affect the accuracy of the flow measurement.

NOTE! DO NOT place the combined CO₂ sensor between the ET tube and the elbow (pediatric/ adult circuit), as this may allow patient secretions to block the adapter windows.

NOTE! The white-striped tubing of the flow sensor should always be in close proximity to the patient.

NOTE! Position the combined CO₂ sensor with its windows in a vertical and NOT a horizontal position: this helps keep patient secretions from pooling on the windows.

NOTE! To prevent rain out and moisture from draining into the tubing of the flow sensor or combined CO₂ sensor, keep the tubing in an upright position.

NOTE! Periodically check the flow sensor and tubing for excessive moisture or secretion buildup.

NOTE! The reading will be a little change when the gas sampling line is mightily occluded. At the same time, this occlusion will damage the module.

Chapter 2: Intended Use and General Information

2.1 Intended Use

The Vital Signs Monitor is a low cost monitor for spot checking, continuous, noninvasive monitoring or recording of functional oxygen saturation of arterial hemoglobin (SpO₂), pulse rate, pulse strength, end-tidal CO₂, FiCO₂ and Respiration Rate. The monitor is a battery or AC adapter powered device. It may be used in the hospital, clinical environment, homecare, and during emergency land transportation.

This device is intended for continuous patient monitoring with adjustable alarm limits as well as visible and audible alarm signals.

NOTE! The monitor was not designed or tested to be an apnea monitor.

2.2 Monitor Features

- Provides fast, reliable SpO₂, pulse rate, pulse strength, ET/Fi CO₂ and RR measurements.
- Ideally suited for use in intensive care units, outpatient clinics, emergency rooms, and during emergency land transport.
- Portable and lightweight. It weighs only 500 grams, with batteries.
- Ergonomically designed to fit comfortably in the palm of your hand.
- Uses four (4) standard “AA” alkaline or Ni-MH batteries.
- Battery life is approximately six (6) hours.
- Bright, easy-to-read LCD displays indicate end-tidal carbon dioxide concentration, SpO₂, pulse rate measurements, Plethysmogram, and trend table.
- Screen rotation provides upright display for vertical or horizontal monitoring positioning.
- Perfusion Index indicates arterial pulse signal strength.

- Two measuring modes:
Spot-Check mode: Measure data intermittently.
Monitor mode: Measure and store data continuously.
- Adjustable volume (including silence) “beep” sounds with each pulse beat.
- Automatic SpO₂, pulse rate, and ETCO₂ alarms with adjustable high and low alarm limits for SpO₂, pulse rate and ETCO₂ measurements.
- Adjustable volume for alarm and alert tones (including silence).
- Low battery icon flashes when about 15 minutes of battery use remains. A red high priority alarm information turns on the alarm bar and an audible two groups of 5-beep burst notifies the user low battery life.
- Patient information management. Patient’s information such as ID, Sex, Type can be managed.
- Data can be transferred to PC through a data line for storage, review and analysis.

2.3 Theory of Operation

2.3.1 Oxygen Saturation (SpO₂)

The device determines the SpO₂ and pulse rate by passing two wavelengths of low intensity light, one red and one infrared, through body tissue to a photodetector. Information about wavelength range can be especially useful to clinicians. Wavelength information for this device can be found in the *SpO₂ Specifications* section of this manual.

Pulse identification is accomplished by using plethysmographic techniques, and oxygen saturation measurements are determined using spectrophotometric oximetry principles. During measurement, the signal strength resulting from each light source depends on the color and thickness of the body tissue, the sensor placement, the intensity of the light sources, and the absorption of the arterial and venous blood (including the time varying effects of the pulse) in the body tissues.




Figure 2.1: Theory of Operation

1. Low intensity Red and Infrared LED light sources

2. Detector

Oximetry processes these signals, separating the time invariant parameters (tissue thickness, skin color, light intensity, and venous blood) from the time variant parameters (arterial volume and SpO₂) to identify the pulses and calculate functional oxygen saturation. Oxygen saturation calculations can be performed because blood saturated with oxygen predictably absorbs less red light than oxygen-depleted blood.

WARNING! Since measurement of SpO₂ depends on a pulsating vascular bed, any condition that restricts blood flow, such as the use of a blood pressure cuff or extremes in systemic vascular resistance, may cause an inability to determine accurate SpO₂ and pulse rate readings.

WARNING! Under certain clinical conditions, the device may display dashes if unable to display SpO₂ and/or pulse rate values. Under these conditions, the device may also display erroneous values. These conditions include, but are not limited to: patient motion, low perfusion, cardiac arrhythmias, high or low pulse rates or a combination of the above conditions. Failure of the clinician to recognize the effects of these conditions on the device readings may result in patient injury.

2.3.2 End-Tidal Carbon Dioxide Concentration (ETCO₂)

The principle is based on the fact that CO₂ molecules absorb infrared (IR) light energy of specific wavelengths, with the amount of energy absorbed being directly related to the CO₂ concentration. When an IR beam is passed through a gas sample containing CO₂, the electronic signal from the photodetector (which measures the remaining light energy) can be obtained. This signal is then compared to the energy of the IR source and calibrated to accurately reflect CO₂ concentration in the sample. To calibrate, the photodetector's response to a known concentration of CO₂ is stored in the monitor's memory. A reference

channel accounts for optical changes in the sensor, allowing the system to remain in calibration without user intervention.

- **2.3.2.1 Mainstream CO2**

The device includes the capability of monitoring end-tidal CO2 using the optional CO2 sensor. This measures CO2 by using the infrared absorption technique, which has endured and evolved in the clinical setting for over two decades and remains the most popular and versatile technique today.

- **2.3.2.2 Sidestream CO2**

The device includes the capability of monitoring end-tidal CO2 using the optional CO2 module, gas dryer line and nasal sampling line.

WARNING! The CO2 module can only supply value of ETCO2, INCO2, respiration rate etc. It is only an assistant to a clinical doctor. The doctor must do the clinical diagnosis according to more information, such as the symptom of the patient, the sickness of the animal, etc.

Chapter 3: Controls and Features

3.1 Monitor Front Panel




Figure 3.1: Monitor Controls, and Features

1. Sensor Connector


The sensor connects here, or an extension cable can be connected between the monitor and the sensor.


2. Pulse Strength Bar Graph


The pulse strength bar graph “sweeps” with the patient’s pulse beat. The height of the bar graph shows the patient’s pulse strength.

3. Mute icon

The mute icon is displayed at the status bar and it has three statuses:

“” this icon means the normal status of alarm sound.

“” this icon is displayed during temporary 30sec, 60sec, 90sec, 120sec alarm silence.

“” this icon is displayed steadily during permanent alarm silence.

4. Power Indicator

This indicator lights steadily to inform the working status of the monitor. Green means the monitor working normally and red means alarm occurred.

5. Pulse Rate Numeric Display

A number shows the patient’s pulse rate value in beats per minute. Dashes (- - -) mean the monitor is not able to calculate the pulse rate value

6. Speaker

It provides audible indication of alarm conditions, pulse tone and feedback for key-presses. Ensure the speaker is not covered.

7. Respiration Rate Numeric Display

Values show the patient's number of breaths per minute. Dashes (---) mean the value of the breathing rate is not available.

8. Information Bar

The information bar displays patient's ID/ type, battery level icon, date/time.

9. Battery level Icon

This icon is displayed at the information bar and has four levels. It flashes when there is only 15 minutes left for the monitor shut down itself.

10. Status Bar

The status bar displays the measuring modes, sensor off/finger off/pulse search/low perfusion icon and volume icon.


11. SpO₂ Numeric Display

A number shows the patient's SpO₂ value in percent. Dashes (- - -) mean the monitor is not able to calculate the SpO₂ value.

12. End-Tidal CO₂ Numeric Display

Value represents the value of the end of the patient's ETCO₂. Dashes (---) mean the system was unable to detect the ETCO₂ value.

End-Tidal CO₂ waveform is shown below:



Phase I : The AB segment baseline (zero), is the beginning of a breath

The phase II : BC segment is the ascending expiratory branch (alveolar and dead space gas mixture)

Phase III: The CD segment breath plateau, horizontal line, mixed alveolar gas

Phase IV: DE segment is the decreased expiratory branch, with a rapid and steep decline to the baseline (fresh gas in the airway)

13. Silence Indicator

This indicator flashes during temporary two-minute alarm silence. The indicator lights steadily during permanent/indefinite alarm silence.

14. USB Interface

The USB interface is used to connect the monitor with PC for trend data output.

15. AC Power Jack

An optional AC power supply connects here.

16. Slot for hanging strap

17. CO2 module state

CO₂✓ Module has been activated.

CO₂✗ Module has not been activated.

3.2 Monitor Operating keys




Figure 3.2: Monitor Operating Keys

1. ON/OFF key

Pressing this key for 5 seconds turns the monitor ON and OFF.

2. Silence key

Pressing the Silence key once in turn can disable the alarm tone for 30sec, 60sec, 90sec, 120sec and disable it indefinitely (until canceled or the monitor is turned off). Note: each pressing should be within 3 seconds.

To cancel the temporary alarm and alert tone silenced condition, press the Silence key twice. To cancel the indefinite silenced alarm, press the Silence key once. The Silence indicator will turn off.

3. Up and Down Arrows

The Up and Down arrow keys are used to adjust the following settings:

- Alarm/ Pulse Volume
- Move the cursor circularly.
- Increase/decrease numbers

- Choose options.

4. Mode Key

Press this key to switch between the four display modes that are big display mode, waveform mode, trend table mode and horizontal display mode.

5. Menu Key

Press this key to change the settings like: patient's information, high/low alarm limits, time and date.

Chapter 4: Operating Instructions

WARNING! Do not use a monitor, sensor, cables, or connectors that appear to be damaged.

WARNING! This device must be used in conjunction with clinical signs and symptoms. This device is only intended to be an adjunct in patient assessment.

4.1 Unpacking the Monitor

1. Carefully remove the monitor and its accessories from the shipping carton. Save the packing materials in case the monitor must be shipped or stored.
2. Compare the packing list with the supplies and equipment you received to make sure you have everything you'll need.

4.2 Install the Batteries




Figure 4.1: Installing the Batteries

The machine uses 4 (four) standard “AA” alkaline or Ni-MH cells.

To install/replace the batteries:

1. Depress the battery door and remove it downward.
2. Install the negative end of each battery first, compressing the battery

terminal spring until the positive terminal clears the positive tab. Press the battery down into place.

3. Place battery door into the slots of the monitor back panel, depress the door tab, and press the door into place.

NOTE! If you install disposable batteries, be sure to dispose of them in compliance with your institution's guidelines and local ordinances.

NOTE! The unit will hold data for about one and a half minutes with no battery power. This will insure the safety of trend data during battery replacement.

4.3 AC Power Adapter

The AC power adapter can be used as the monitor's power supplier and Ni-MH cells' charger.

NOTE! Do not plug the monitor into an outlet controlled by a wall switch.

NOTE! When using AC power, the digital device is a class II device with functional earth. This earth connection is for device electromagnetic compatibility and does not provide protection to the patient or user.

NOTE! Do not use the AC adapter to charge alkaline batteries.

4.3.1 Charging Ni-MH batteries

The battery may discharge during prolonged storage or shipment. If the monitor have been in storage for more than 2 months, it is important to plug the AC power adapter into an AC outlet and allow the batteries to charge for approximately 30 minutes before attempting to operate the instrument.

To charge a low battery, connect the monitor to an AC power through the AC power adapter. A full charge of a completely discharged battery takes 20 hours while the monitor is turned off.

4.4 Attaching the Sensor to the Patient

WARNING! Prolonged use or the patient's condition may require changing the sensor site periodically. Change sensor site and check skin

integrity, circulatory status, and correct alignment at least every 4 hours.

Attaching the patient to the monitor requires these steps:

1. Choose SpO₂ and ETCO₂ sensors.
2. Check the sensors and oximetry cable.
3. Clean or disinfect the sensor if using the reusable type (Disposable sensors are for single-patient use and do not require cleaning or disinfecting).
4. Attach the sensor to the patient.

4.4.1 Choosing the Sensor

WARNING! Before use, carefully read the sensor directions for use, including all warnings, cautions, and instructions.

Choose the appropriate SpO₂ sensor from the following chart.

PATIENT	SITE	DESCRIPTION
Adult >45Kg	Finger	Sensor, Adult (reusable)
	Finger or toe	Sensor, Disposable, Adult Finger
	Ear	Sensor, Ear (reusable)
Pediatric 15-45Kg	Finger	Sensor, Adult (reusable)
	Finger or toe	Sensor, Disposable, Ped. Finger
	Ear	Sensor, Ear (reusable)
Neonate <3Kg (for spot-check only)	Hand or Foot	Sensor, Disposable, Neonate
	Foot	Sensor, Wrap. Neonate(reusable)

4.4.2 Care and Handling of the Sensor

WARNING! Misuse or improper handling of the sensor and cable could result in damage to the sensor. This may cause inaccurate readings.

Hold the connector rather than the cable when connecting or disconnecting the sensor to the device as shown in Figure 4.2.

Adapter Cable




Figure 4.2: Disconnecting or connect the sensor

There are three connectors on the adapter cable. The one is used to plug on the device. The other two connectors are used to connect the SpO₂ and ETCO₂ sensor.

NOTE! Please plug the SpO₂ and ETCO₂ sensors on its own connectors correctly. Otherwise, the device may work abnormal.

Do not use excessive force or unnecessary twisting when connecting, disconnecting, storing, or when using the sensor.

Placing the adult/pediatric SpO₂ sensor:

When placing the sensor on the patient, allow the cable to lay the back of hand as shown in Figure 4.3.




Figure 4.3: Positioning the cable of the finger sensor

Placing the neonatal SpO₂ sensor:

Neonatal SpO₂ sensors are Y-type with a rubber wrap. First insert the Y into the slots of the wrap as shown in Figure 4-4. After placed, the Neonatal SpO₂ sensors look like Figure 4-5.




Figure 4.4 Placing Neonatal SpO₂ Sensor




Figure 4.5 Placing Neonatal SpO₂ Sensor

Place the sensor onto the foot or hand. Secure the wrap (about 20mm long) to ensure the correct position of the sensor as shown in Figure 4-6. Do not secure the wrap too tight as it may affect the blood flow.




Figure 4.6 Placing Neonatal SpO₂ Sensor

Placing the Mainstream ETCO₂ Sensor:

The CO₂ sensor is a rugged, solid-state, mainstream sensor. It is factory calibrated and does not require further calibration.




Figure 4.7 ETCO₂ Sensor

Insert the CO₂ sensor connector into the adapter cable's connector which was

printed “CO2 Sensor”.

To remove the connector, grasp the body portion of connector back and remove.

Shown below is the CO2 sensor with a patient circuit:




Figure 4.8 Placing Neonatal SpO₂ Sensor

Placing the Sidestream ETCO₂ sensor:

The CO2 module is a rugged, solid-state, sidestream sampling mode. It is factory calibrated and does not require further calibration.



Figure 4.7: CO2 module

- 1). Insert the CO2 module connector into the adapter cable’s connector which was printed “CO2 Sensor”.
- 2). Insert the gas dryer line to CO2 module first, then insert the nasal sampling line to gas dryer line.
- 3). Attaching the nasal sampling line to patient, shown as figure 4.8:




Figure 4.8: Placing the Nasal Sampling Line.

WARNING! The nasal sample line can be used just once. If you use it repeatedly, it may cause cross infection.

CAUTION! The Gas Dryer line can be used only 150 hours. If use it too long, it may cause wrong data or damage the module.

CAUTION! Do not use the module without connecting with the Drying line. Otherwise, it will cause measurement error or damage the module.

CAUTION! The measurement will be error when the module is not warm up.

NOTE! The reading will be not accurate when using in hard wavy-temperature environment.

NOTE! The reading will be a little change when the gas sampling line is mightily occluded. At the same time, this occlusion will damage the module.

NOTE! Gas leakage by any reason will cause the reading fall down and influence the configuration of CO₂ wave.

NOTE! Please use the accessories supplied by the manufacturer, otherwise, the reading may not be accurate.

4.4.3 Checking the Sensor and Adapter Cable

Follow these instructions each time before you attach the sensors to the patient. This helps ensure the sensors and adapter cable are working properly.

WARNING! Using a damaged sensor/cable may cause inaccurate readings. Inspect each sensor/cable. If a sensor/cable appears damaged, do not use it. Use another sensor/cable or contact your authorized repair center for help.

1. Carefully inspect the SpO₂ and ETCO₂ sensor to make sure it does not appear damaged.
2. Carefully inspect the adapter cable to make sure it does not appear damaged.
4. When using the adapter cable:
 - a. If the sensors are not already connected to the cable, connect the sensors to the cable firstly. Push the connectors together firmly and close the latch to secure the connectors.

- b. If the cable is not already connected to the monitor, connect the adapter cable to the monitor. Push the connector firmly into the monitor.
5. If not using the adapter cable, connect one of the SpO₂/ETCO₂ sensor to the monitor. Push the connector firmly into the monitor.
6. If the monitor is not already on, press the On/Off key to turn on the monitor.

WARNING! If any of the integrity checks fail, do not attempt to monitor the patient. Use another sensor or oximetry cable, or contact the equipment dealer for help if necessary.

7. Before the sensor is attached to the patient, check the integrity of the sensor, adapter cable.
 - a. Make sure the red light in the sensor is illuminated.

NOTE! Obstructions or dirt on the sensor's red light or detector may cause the checks to fail. Make sure there are no obstructions and the sensor is clean.

NOTE! You can NOT measure the ETCO₂ in 2 to 3 minutes after the monitor turning on. Because the monitor is in the process of zeroing during this period.

8. The ETCO₂ module will conduct zeroing automatically after warming up for 2 to 3 minutes.

NOTE! Please make sure there is no CO₂ gas existing in the airway adapter and conduct the ZEROING process button before using. Otherwise, the measurement may be inaccuracy.

NOTE! The failure of zeroing may be caused by:

- In sleep mode
- CO₂ module temperature not stable
- CO₂ module sensor faulty
- CO₂ module Zero attempted and breaths have been detected in the last 20 seconds.

9. After the warming up and zeroing, you are now ready to attach the sensors to the patient.

4.4.4 Cleaning or Disinfecting the Sensors

Clean or disinfect reusable sensors before attaching to a new patient.

WARNING! Do not autoclave, ethylene oxide sterilize, or immerse the sensors in liquid.

CAUTION! Unplug the sensor from the monitor before cleaning or disinfecting.

Clean the sensor with a soft cloth moistened in water or a mild soap solution. To disinfect the sensor, wipe the sensor with isopropyl alcohol.

4.5 Performance Considerations

4.5.1 SpO₂ Performance Considerations

WARNING! SpO₂ readings and pulse signal can be affected by certain ambient environmental conditions, sensor application errors, and certain patient conditions.

Inaccurate measurements can be caused by:

- Incorrect application of the sensor
- Placement of the sensor on an extremity with a blood pressure cuff, arterial catheter, or intravascular line
- Ambient light
- Prolonged patient movement

Loss-of-pulse signal can occur for the following reasons:

- The sensor is too tight
- A blood pressure cuff is inflated on the same extremity as the one with the sensor attached
- There is arterial occlusion proximal to the sensor

Select an appropriate sensor, apply it as directed, and observe all warnings and Cautions presented in the directions for use accompanying the sensor. Clean and remove any substance such as nail polish from the application site. Periodically check to ensure that the sensor remains properly positioned on the patient.

WARNING! Tissue damage can be caused by incorrect application or duration of use of a SpO₂ sensor. Inspect the sensor site as directed in the sensor directions for use.

High ambient light source such as surgical lights (especially those with a xenon light source), bilirubin lamps, fluorescent lights, infrared heating lamps, and direct sunlight can interfere with the performance of a SpO₂ sensor. To prevent interference from ambient light, ensure that the sensor is properly applied, and cover the sensor site with opaque material.

NOTE! Failure to take this precaution in high ambient light conditions may result in inaccurate measurements.

If patient's movement presents a problem, try one or more of the following remedies to correct the problem:

- Verify that the sensor is properly and securely applied
- Move the sensor to a less active site
- Use an adhesive sensor that tolerates some patient motion
- Use a sensor with fresh adhesive backing

4.5.2 ETCO₂ Performance Considerations

- a. Inaccurate measurements for mainstream ETCO₂ can be caused by:
 - Forget to warm up the monitor in 2 to 3 minutes.
 - CO₂ exists in the airway adapter. Please make sure there is no CO₂ exist in the airway adapter and conduct the ZEROING process button before using.
 - Forget to set the CO₂ settings. You need first adjust all the ETCO₂ settings correctly before using.
 - The CO₂ sensor is too far away from the patient. Nearer is better.
- b. Inaccurate measurements for sidestream ETCO₂ can be caused by:
 - Forget to set the CO₂ settings. You need first adjust all the ETCO₂ settings correctly before using.
 - Working at strong disturbing of electromagnetism (such as ESU, MRI, CT etc.) environment, the module will be inaccurate.
 - Using it in hard wavy-temperature environment, the module will be inaccurate.
 - The reading will be a little change when the nasal sampling line is mightily occluded. At the same time, this occlusion will damage the module.
 - Gas leakage by any reason will cause the reading fall down and influence the configuration of CO₂ wave.


4.6 Turning On the Monitor

WARNING! To ensure patient safety, do not place the monitor in any position that might cause it to fall on the patient.

WARNING! As with all medical equipment, carefully route patient cabling to reduce the possibility of patient entanglement or strangulation.

WARNING! Ensure that the speaker is clear of any obstruction. Failure to do so could result in an inaudible alarm tone.

NOTE! You can NOT measure the ETCO₂ in 2 to 3 minutes after the monitor turning on. Because the monitor is in the process of zeroing and the ETCO₂ sensor needs to warm up.

1. To turn on the monitor, press and held the ON/OFF key for about five seconds. When turned ON, the monitor does the following:
 - The monitor starts a power-on self test (POST) quickly.
 - The monitor's software revision is momentarily displayed.
 - The Power Indicator illuminates green.
2. Upon successful completion of the POST, the monitor sounds a one-second tone indicating that the monitor has passed the test.
3. If a sensor is connected to the monitor and the patient, the Pulse Search Icon  is displayed at the Status Bar. "--" will flash on the numeric display until the SpO₂ and pulse rate reading have satbilized. The monitor will search the pulse for 8 seconds.
4. The monitor will pop up a dialog to remind operator to set the patient's ID. If choose NO or no pressing within 5 seconds, the monitor will use the last ID for the present patient.
5. Monitor the patient.


WARNING! Verify that the power indicator lights up and you can hear the POST pass tone upon startup of the device. If not, do not use the monitor.

WARNING! The monitor will automatically be powered off when no finger is in the device and no operation for longer than five minutes in the Spot-check and Monitoring measuring modes. The screen brightness will be decreased when no finger is in the


device and no operation for longer than three minutes in the Recording measuring modes.

4.6.1 Six Display Modes


There are six Display Modes, you can switch between them by pressing the Mode Key. The modes are shown as below:




(A) Big Display Mode




(B) Waveform Display Mode




(C) Trend-table Display Mode



(D) Horizontal Display Mode



(E) SpO2 Figure Mode



(F) ETCO2 Trend Diagram Mode

Figure 4.9 Six Display Modes

4.6.2 Two Measuring Modes

There are two Measuring Modes: Spot-Check and Continuous Monitor. Their differences are compared in to below table:

Modes	Spot-Check	Monitor
Way of Measurement	intermittently	continuously
Data Saved or not	NO	YES
Energy-save or not	NO	NO
Alarm or not	YES	YES
Volume Adjustable or not	YES	YES

4.7 Turning Off the Monitor

After monitoring is complete, remove the sensor from the patient and store or dispose of the sensor according to governing rules.

Turn off the monitor when you are not monitoring a patient. To turn off the monitor, press and held the On/Off key for about five seconds.

Chapter 5: Changing the Monitor's Settings

5.1 Changing the Pulse Volume



Figure 5.1 Pulse Volume Bar

A “beep” tone sounds with each pulse beat. The volume of the “beep” can be adjusted to seven (7) settings. Setting the volume by pressing the Up and Down Arrow, the volume is increased or decreased with each key press.

5.2 Changing the Alarm Volume



Figure 5.2 Alarm Volume Bar

Hold the Silence key for 3 seconds to change the alarm volume. Setting the volume by pressing the Up and Down Arrow, the volume is increased or decreased with each key press.

5.3 Managing the Patient's Information

Please write down the information of the patient who is going to be monitored like: Patient's ID (from 000 to 999), Sex (Male or Female), Type (Adult/ pediatric /neonatal) and Mode (Spot-Check/Record/Monitor).

Main Menu		
PATIENT	ID	678
SPO2 ALM	Sex	Male
	Type	Adu
CO2 ALM	Mode	Spot-Check
CO2 SET		
SYS SET		
RESET		
		OK Cancel




Figure 5.3 Managing the patient's information

1. Press the Menu key to enter into the Main Menu interface.
2. Press the Down key and the Menu key to enter the "Patient" submenu.
3. Press the Up/Down key to choose items you want to change and press the Menu key to active the items.
4. Press the Up/Down key increase/decrease the setting and then use the Menu key again to exit.
5. Press the Up/Down key to OK/ Cancel button to confirm/ cancel your settings.

5.4 Changing the Alarm Limits

Alarms are audio and visual indicators generated by the monitor to alert doctors and nurses. These alarms occur when the vital signs of the patients being monitored become abnormal, or the monitor itself malfunctions and could not perform the monitoring task.

5.4.1 Alarm Summary

Alarm Events	Alarm Level	Alarm adjustable or not	Visual alarm	Audio alarm
SpO ₂ exceed the pre-set limits	High	Yes	<ol style="list-style-type: none"> 1. Numbers of SpO₂ will flash. 2. Alarm event will be displayed in red at the Alarm Bar. 3. Alarm icon  will flash in the middle of the screen for 5 seconds. 	Two groups of "dee-dee-dee-dee-dee" Period: 10 secs
PR exceed the pre-set limits	Medium	Yes	<ol style="list-style-type: none"> 1. Numbers of SpO₂ will flash. 2. Alarm event will be displayed in yellow at the Alarm Bar. 	Two groups of "dee-dee-dee" Period: 18 secs
Sensor off	High	NO	<ol style="list-style-type: none"> 1. The icon  will be displayed at the Status Bar. 2. "Sensor off" will be displayed in red at the Alarm Bar. 	Two groups of "dee-dee-dee-dee-dee" Period: 10 secs
Finger off	High	NO	<ol style="list-style-type: none"> 1. The icon  will be displayed at the Status Bar. 2. "Finger off" will be displayed in red at the 	Two groups of "dee-dee-dee-dee-dee" Period: 10 secs

			Alarm Bar.	
Low battery	High	NO	1. The Battery Level Icon will flash. 2. "Low battery" will be displayed in red at the Alarm Bar.	Two groups of "dee-dee-dee-dee-dee" Period: 10 secs
Data full	Medium	NO	"Data full" will be displayed in yellow at the Alarm Bar.	Two groups of "dee-dee-dee" Period: 18 secs
Battery abnormal	High	Yes	"Battery abnormal" will be displayed in red at the Alarm Bar.	Two groups of "dee-dee-dee-dee-dee" Period: 10 secs
ETCO2 exceeds the pre-set limits	Medium	Yes	1. ETCO2 value flashes 2. Alarm event will be displayed in red at the Alarm Bar. 3. Alarm event will be displayed in yellow at the Alarm Bar.	Two groups of "dee-dee-dee" Period: 18 secs
RR exceeds the pre-set limits	Medium	Yes	1. RR value flashes 2. Alarm event will be displayed in red at the Alarm Bar. 3. Alarm event will be displayed in yellow at the Alarm Bar.	Two groups of "dee-dee-dee" Period: 18 secs

NOTE! The power indicator will flash in red in all levels of alarm.

WARNING! There is no alarm in the Record Measuring Mode.

5.4.2 Changing the SpO2 Alarm Settings




Figure 5.4 Changing the SpO2 alarm settings

Each measurement, SpO₂ and Pulse Rate, has a high and low alarm limit setting.

1. Press the Menu key to enter into the Main Menu interface.
2. Press the Down key twice and the Menu key to enter the “SPO2 Alarm” submenu.
3. Press the Up/Down key to choose items you want to change and press the Menu key to start changing the settings of a particular parameter.
4. Press the Up/Down key increase/decrease the setting and then use the Menu key again to exit.
5. Press the Up/Down key to the OK/ Cancel button to confirm/ cancel your settings.

You can press the “Reset” button to reset you settings.

5.4.3 Changing the ETCO2 Alarm Settings




Figure 5.5 Changing ETCO2 and respiration rate alarm settings

ETCO₂ and respiratory rate measurements also need to set upper and lower alarm limits.

1. Press the Menu key to enter into the Main Menu interface.
2. Press the Down key three times and the Menu key to enter the "CO₂ Alarm" submenu.
3. Press the up/down keys to select the appropriate item, and then press the Menu key to enter the setup.
4. Press up/down keys settings to modify settings, and then press the Menu key to exit.
5. Press up/down keys to save or cancel the setup.

You can press the reset button to reset your settings.

WARNING! Be aware of alarm limits of similar units in the same area when adjusting alarm limits of this device to avoid confusion.

NOTE! Alarm limits are non-overlapping. You cannot set the high alarm equal to or lower than the low alarm and you cannot set the low alarm equal to or higher than the high alarm.

NOTE! While setting alarm limits, if no keys are pressed for twenty seconds, the alarm limit setting mode is exited and the SpO₂ and pulse rate measurements are shown. Changes aren't saved.

NOTE! The alarm actions occur for each violated alarm, even if more than one alarm is violated at the same time.

NOTE! Alarms may be tested while the monitor is in use by setting alarm limits such that the measured parameter is outside alarm limits. Return limits to the required settings after testing.

NOTE! If you change the patient's type, the alarm limits will become default limits of this type.

5.5 Changing ETCO₂ Settings

CO₂ Compensations for Mainstream and Sidestream

The measurement of CO₂ is affected by temperature, pressure, and gas compensations. The barometric pressure is as well as the presence of O₂,

N₂O, helium, and anesthetic agents in the gas mixture need to be compensated for by the device in order to achieve its stated accuracy. The device provides instrument settings to allow the host to communicate these operating conditions. The instrument settings for these parameters should be set when initially connecting to the device and whenever there is a change in the conditions at the patient airway.

The following table shows the affect of various concentrations of expired gas on the reported CO₂ values if the default compensations are not adjusted to match the mixture. The table assumes a CO₂ value of 38 mmHg, a default O₂ setting of 16%, all other compensations are 0%. The airway gas is 5% CO₂, 16% O₂, the temperature is 35 °C, and the barometric pressure is 760 mmHg.

Table: Gas Compensation Effects on CO₂

Expired Gas Mixture				
% O ₂	Balance Gas	Baro Pressure mmHg	Gas Temp °C	Reported CO ₂ mmHg
16	N ₂	760	35	38
60	N ₂	760	35	36.6
60	N ₂ O	760	35	39.9
40	N ₂ O	760	35	40.3
16	N ₂	700	35	33.9 ⁽¹⁾
16	N ₂	760	25	39.1

Notes:

(1) At 700mmHg of pressure, the correct CO₂ value is 35.0 mmHg




Figure 5.6 Changing ETCO₂ Settings

Please set the correct CO2 settings according to your operation environment at the first time you using this monitor.

1. Press the Menu key to enter into the Main Menu interface.
2. Press the Down key three times and the Menu key to enter the "CO2 SET" submenu.
3. Press the up/down keys to select the appropriate item, and then press the Menu key to enter the setup.
4. Press up/down keys settings to modify settings, and then press the Menu key to exit.
5. Press up/down keys to save or cancel the setup.

You can press the reset button to reset your settings.

5.6 System Settings

5.6.1 Changing the ETCO2 Alarm Settings




Figure 5.8 Changing the date and time

1. Press the Menu key to enter into the Main Menu interface.
2. Press the Down key thrice and the Menu key to enter the "SYS SET" submenu.
3. Press the up/down keys to select the appropriate item, and then press the Menu key to enter the setup.
4. Press up/down keys settings to modify settings, and then press the Menu key to exit.
5. Press up/down keys to save or cancel the setup.

You can press the reset button to reset your settings.

To prolong the life of CO2 module, we designed the option of "CO2 Active".

There are 3 options in "CO2 Active": Auto, Enable, Disable.

- 1) Auto: CO2 can switch to working mode and sleep mode automatically. When the SpO2 and CO2 don't measurement about 3 minutes, CO2 module will switch to sleep mode automatically.

NOTE! It is unable to measurement CO2 in sleep mode. Until the SpO2 be measured or any other operates, then CO2 will switch to working mode.

- 1) Enable: CO2 is in work mode, and it is measuring.
- 2) Disable: CO2 is in sleep mode, unable to measurement.

5.7 Resetting




Figure 5.9 Restoring default settings

1. Press the Menu key to enter into the Main Menu interface.
2. Press the Down key for four times and the Menu key to select the "Reset" option.

NOTE! If you choose to reset your settings, all settings will be restored to the Factory Default Settings except the Date and Time.

Factory Default Settings:				SpO ₂	PR	ETCO ₂	RR	
Patient	ID	987	Alarm	Alarm	On	On	On	On
	Sex	Male		Upper	99	100	40	30
	Type	Adult		Lower	92	60	30	8
	Mode	Spot- Check		Level	H	M	M	M

Chapter 6: Trend Data Output and Analysis

6.1 Connecting the Device to PC

The device can store 4750 hours of SpO₂, Pulse Rate, EtCO₂, FiCO₂, RR, Trend Data captured at 1 minute intervals. This trend data can be transferred to a PC for evaluation.

Trend data is stored in non-volatile memory, so it is not erased when the unit is shut off or when the batteries are replaced.

The device will send out trend data through USB interface so that the data can be stored, analyzed and printed. An additional USB Data Cable is required to connect the device to a PC.




Figure 6.1: PC Communication Setup

6.2 Installing the Data Management System

The Data Management System (DMS) is a management and analysis system for history data. It supplies functions like data download, data review, case information edit, data analysis and result print etc. To install the DMS, Please conduct the following steps:

1. Insert the CD into CD-ROM.
2. Copy to the Setup file from the CD to your computer.
3. Open the Setup file and doubleclick  to install it. It will take about 1 minute.
4. You are now ready to use the DMS.

6.3 Function of Main Menu




Figure 6-2: The main menu



1. **Open Data Base:** Allows user to create, delete, search, and view patient files and study records.



2. **Case Information:** Allows user to immediately refer back to patient information while viewing study data.



3. **Down Load All:** Allows user to access to the patient database and transfer all trend data from the device to ODMS for analysis.



4. **Print Report:** Allows user to print statistical, desaturation, SpO₂, or full reports.



5. **Option:** Allows user to setup language etc.



6. **Help:** Displays user sheet and operate method when problems occur.




7. **Exit:** Allows user to exit the DMS.


6.4 Trend Data Download

To download trend data to PC, please conduct the following steps:

1. Connect the device with the PC through an USB data cable.
2. Turn on the device.

3. Start the DMS.
4. Select the appropriate COM port number, if necessary.
5. There are two way to download trend data:

A. Click the icon “” or “Case Info” at the main menu, the following interface will pop up:



The screenshot shows a window titled "Current Case Information" with a close button in the top right corner. The window is divided into several sections:

- Patient:** Contains three text input fields: "First Name" (value: Peter, <33 bytes), "Last Name" (value: Jiang, <33 bytes), and "Medical Record ID" (value: 985, xxx(001)). Below these is a "Detail>" link and a folder icon.
- Case Date/Time:** Contains two text input fields: "Date:(mm/dd/yy)" (value: 2010-01-26) and "Time:(hh:mm:ss)" (value: 18:10:26).
- DownLoad Date/Time:** A text field showing the value 2010-01-29 12:02:49.
- Medical Information:** Contains two text areas: "Indication" and "Comments".
- Physician Name:** A text input field.
- Buttons:** At the bottom right, there are three buttons: "Begin DownLoad" (with a small icon), "Save", and "Cancel".


Figure 6.3: Case Information

Patient: Fill in the basic and detail information of the patient.


Case Date/Time: Means the time to collect the patient’s data.


DownLoad Date/Time: Means the time to stop collecting the patient’s data.

Medical Information: Fill in the name of physician, indication of the patient and comments on the patient.


After writing down the current case information, click this  button to begin download the trend data of this ID.

NOTE! The patient’s ID and information should be corresponded with each other.

If you want to change the case information of the patient whose trend data have already been download, follow the same steps and click this  button to save your changes.

B. Click the icon “” or “Download All” to download all trend data from the device.

6.5 Open the Data Base

Click the icon “” or “Open DataBase” to open the data base, the following window will pop up:






Figure 6.4: Data Base

Case Search: Put ID/First Name/Last Name to search patient files you want to review.

Case List: List all the downloaded cases. “???” means the information of this patient hasn’t been edited. Click the  button to add the information.

Case Test List: Double click the patient file, the case test list of this patient will be displayed at the right side.

Delete: Select patient file you want to delete and then click the  button to delete it.

6.6 Data Analysis and Reports

Select one period from the case test list to analyze trend data. Press the OK button to enter the main analysis interface, as shown below:




Figure 6.5: The Main Analysis Interface

The main analysis interface displays two bar charts and one list of SpO₂ and PR.

The Full Report consists of 3 reports: Statistical Report, Desaturation Report and Graphic Report.

1. **Statistical Report:** The software will analyze, calculate, and print statistical parameters for downloaded data.

Chapter 6: Trend Data Output and Analysis




Figure 6.6: The Statistical Report

Screening Start Date/Time: represents when the monitor began collecting data.

Screening End Date/Time: Time when monitor stopped collecting data, either because of power off or memory full, and is calculated by the software based on the amount of data downloaded.

Total Screening Duration: calculated by software.

Total Valid Sample Duration: Total Screening duration minus the excluded sample time.

Total Excluded Sample Time: Total time the monitor registered invalid data.

Low SpO₂: The lowest SpO₂ value taken from the valid sample with the pulse rate and time at which the lowest SpO₂ value occurred.

SpO₂ Standard Deviation: The measure of the variability of the sample. A higher standard deviation indicates greater variability in the samples.

High PR: The highest pulse rate value taken from the valid sample and the SpO₂

value & time the highest pulse rate occurred.

Low PR: The lowest pulse rate taken from the valid sample and the SpO₂ value & time the lowest pulse rate occurred.

2. Desaturation Report: Lists up to 50 of the deepest (lowest SpO₂ values) desaturation events. The events are listed in order of depth, beginning with the deepest event.

Sleep Screening Desaturation Report							
The following lists the 48 desaturation events having the lowest SpO ₂ values. They are listed in order by depth of desaturation, beginning with the deepest event. An event is defined as a desaturation greater than or equal to 4% with resaturation greater than or equal to 3%, OR a desaturation greater than or equal to 3% with a resaturation greater than or equal to 4%.							
n/n	Begin Time	Stop Time	Duration (min:sec)	Min SpO ₂ (%)	Average PR(BPM)	PR Max (BPM)	PR Min (BPM)
1	15:35:32	15:35:47	00:00:15	47	36	38	32
2	10:57:49	10:58:08	00:00:20	50	28	28	28
3	12:11:35	12:12:07	00:00:32	65	211	245	183
4	09:41:56	09:42:26	00:00:30	70	116	146	104
5	11:01:20	11:01:39	00:00:19	75	178	187	167
6	16:27:39	16:28:16	00:00:37	76	29	30	29
7	17:51:40	17:51:53	00:00:05	77	29	29	29
8	15:07:06	15:07:12	00:00:06	78	28	28	28
9	17:54:53	17:55:20	00:00:27	78	181	181	181
10	11:46:57	11:47:28	00:00:31	79	146	146	146

Desaturation Summary	
Total number of desaturation events:	48
The average minimum for SpO ₂ desaturation events:	81%
The following lists the periods of time during which the patient's SpO ₂ remained at or below 88% and 85%.	
There were 0 periods with an SpO ₂ 88% or lower, that were 14 minutes in duration or greater.	
The longest single period with an SpO ₂ 88% or below was:	2 minutes and 17 seconds in duration at 16:17:53 time.
There were 0 periods with an SpO ₂ 85% or lower, that were 12 minutes in duration or greater.	
The longest single period with an SpO ₂ 85% or below was:	2 minutes and 13 seconds in duration at 10:08:32 time.

Figure 6.7: The Desaturation Report

Reference for the event: A unique number between 1-50, which identifies each desaturation event.

Duration of the event: calculated by the software.

Desaturation Summary:

- Total number of desaturation events detected
- Average minimum for SpO₂ desaturation events
- Number of events with 88% SpO₂ or below 5 minutes or longer
- Duration of longest single event with 88% SpO₂ or below & time it occurred
- Number of events with 85% SpO₂ or below 5 minutes or longer

- Duration of longest single event with 85% SpO₂ or below & time it occurred

3. **Graphic Report:** it contains three trend graphs of SpO₂,PR,EtCO₂,RR per page.




Figure 6.8: The Graphic Report

Chapter 7: Accessories

7.1 Standard Configuration

a. Mainstream

#	DESCRIPTION	QTY
1	Vital Signs Monitor	Each
2	Adult Finger SpO ₂ Sensor	Each
3	ETCO ₂ Module	Each
4	Adult Airway Adapter	Each
5	Pediatric Airway Adapter	Each
6	DB9 Adapter Cable	Each
7	Hanging Strap	Each
8	Operation Manual	Each

b. Sidestream

#	DESCRIPTION	QTY
1	Vital Signs Monitor	Each
2	Adult Finger SpO ₂ Sensor	Each
3	DB9 Adapter Cable	Each
4	ETCO ₂ Module	Each
5	Gas Dryer Line	Each
6	Nasal Sampling Line	Each
7	On-airway Sampling Line	Each
8	T-fitting	Each
9	Hanging Strap	Each
10	Operation Manual	Each

7.2 Optional Accessories

#	DESCRIPTION	QTY
1	AC Adapter	Each
2	Data Analysis Software CD	Each
3	AA size Ni-MH Batteries	Four
4	Reusable sensor, Adult, Finger	Each
5	Reusable sensor, Pediatric, Finger	Each
6	Reusable sensor, Neonatal/infant, Foot	Each
7	Single patient use sensor, Adult>30kg	Each
8	Single patient use sensor, Ped 3-50kg	Each
9	Single patient use sensor, neonatal <3kg	Each
10	SpO ₂ extension cable	Each

Chapter 8: Maintenance and Troubleshooting

8.1 Schedule of Maintenance

MAINTAIN THIS ITEM	HOW OFTEN	BY DOING THIS
Battery	When Battery Level icon is flashing, and/or audible alarm sounds.	Follow the instructions for installing the batteries.
Disinfecting the reusable sensor.	Before attaching the sensor to the patient.	Follow the instructions for cleaning the reusable sensor.
Disinfecting the monitor.	When necessary.	<ol style="list-style-type: none"> 1. Remove the batteries from the unit. 2. Wipe the surfaces of the monitor with a soft, clean cloth dampened in isopropyl alcohol. Use only a cloth that is dampened, not wet.

CAUTION! Do not allow isopropyl alcohol or water to enter any of the openings on the monitor. Evidence that liquid has been allowed to enter the monitor voids the warranty.

8.2 Storage

WARNING! To ensure accurate performance and prevent device failure, do not subject the device to extreme moisture such as direct exposure to rain. Such exposure may cause inaccurate performance or device failure.

Whenever possible, the monitor should be stored at room temperature in a dry environment.

If it is necessary to store the monitor for an extended period of time, the unit should be packed in its original shipping container. Storing the monitor for a long period of time may degrade the battery capacity. Batteries should be removed from the monitor before storing.

Storage specifications are as follows:

Temperature: -20 °C to +55 °C

Relative Humidity: 10% to 95% (noncondensing)

8.3 Troubleshooting

PROBLEM	POSSIBLE CAUSE	CORRECTIVE ACTION
No pulse shown on the bargraph.	<ul style="list-style-type: none"> • Patient cable or sensor is disconnected from the device. • Sensor is incorrectly positioned on the patient. • Poor patient perfusion. • Defective sensor or patient cable. 	<ul style="list-style-type: none"> • Check sensor connections to the patient cable and to the device. • Reposition the sensor. • Reposition the sensor. • Try a new sensor or contact your authorized repair center for help.
Pulse rate is erratic, intermittent, or incorrect.	<ul style="list-style-type: none"> • Sensor incorrectly positioned. • Patient motion 	<ul style="list-style-type: none"> • Reposition the sensor. • Patient must remain still to obtain an accurate measurement.
SpO2 value is erratic, intermittent, or incorrect.	<ul style="list-style-type: none"> • Poor patient perfusion. • Patient motion. 	<ul style="list-style-type: none"> • Reposition the sensor. • Patient must remain still to obtain an accurate measurement.
No PR and SpO2 values.	<ul style="list-style-type: none"> • Defective sensor or patient cable or monitor. 	<ul style="list-style-type: none"> • Try a new sensor or contact your authorized repair center for help.
Battery Abnormal	<ul style="list-style-type: none"> • Charging alkaline batteries. • Batteries incorrectly installed. • There are no batteries. 	<ul style="list-style-type: none"> • Install Ni-MH cells instead of alkaline ones. • Reposition batteries correctly. • Equip the device with batteries.
The device doesn't turn on.	<ul style="list-style-type: none"> • Batteries weak. • Batteries not installed or batteries incorrectly installed. 	<ul style="list-style-type: none"> • Replace the batteries. • Ensure the batteries are installed correctly.
The device turns off unexpectedly.	<ul style="list-style-type: none"> • Batteries are weak or dead. 	<ul style="list-style-type: none"> • Replace the batteries.

Sensor	<ul style="list-style-type: none"> •Patient cable or sensor is disconnected from the device. •Sensor is incorrectly positioned on the patient. •Poor patient perfusion. •Defective sensor or patient cable 	<ul style="list-style-type: none"> •Check sensor connections to the patient cable and to the device. •Reposition the sensor. •Reposition the sensor. •Try a new sensor or contact The manufacturer After Service Department for help.
Carbon dioxide waveform disappears, and air pressure decreases	<ul style="list-style-type: none"> • Connector is leaking 	<ul style="list-style-type: none"> • Check the endotracheal tube threaded pipe • Check threaded pipe anesthesia machine or breathing bag •If there is surgery involving the face, make sure the respiratory face mask is secured properly.
Carbon dioxide waveform disappears, and air pressure increases	<ul style="list-style-type: none"> • Catheter is coiled • Airway obstruction • Valve failure 	<ul style="list-style-type: none"> • Remove the catheter obstruction • Remove the airway obstruction • Replace the valve
Sensor Over Temp	<ul style="list-style-type: none"> • The sensor temperature is greater than 50 °C. 	<ul style="list-style-type: none"> • Make sure sensor is not exposed to extreme heat (heat lamp, etc.). If error persists, return sensor to factory for servicing.
Sensor Faulty	<ul style="list-style-type: none"> • ETCO2 Module Source Current Failure • EEPROM Checksum Faulty • Hardware Error 	<ul style="list-style-type: none"> • Check that the sensor is properly plugged in. Reinsert or reset the sensor if necessary. If error persists, return sensor to
In Sleep Mode	<ul style="list-style-type: none"> • This bit is set when sensor has been placed in sleep mode. 	/

Zero In Progress	• Zero is currently in progress.	/
Sensor Warm Up	If calibration is required while this error is set, the message “Wait for Sensor” can be shown.	/
Check Sampling Line	/	Check that the sampling line is not occluded or kinked.
Zero Required	Require the Zeroing process.	/
CO2 Out of Range	Remind message	/
Check Airway Adapter	Remind message	/

8.4 EMC (Electro-Magnetic Compatibility)

Caution! Vital signs monitor complies with the requirements for medical devices EN60601-1-2: 2001, Medical Device Directive 93/42/EEC. Vital signs monitor has passed tests for CISPR 11 class A.

Guidance and manufacturer’s declaration-electromagnetic emissions For all EQUIPMENT and SYSTEMS

Guidance and manufacturer’s declaration-electromagnetic emission		
<i>Vital signs monitor</i> is intended for use in the electromagnetic environment specified below. The user of <i>Vital signs monitor</i> should make sure that it is used in such environment.		
Emission test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	<i>Vital signs monitor</i> uses RF energy only for its internal function. Therefore its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class A	<i>Vital signs monitor</i> is suitable for use in all establishments other than domestic and those
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	applicable	Connect to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Guidance and manufacturer's declaration-electromagnetic immunity For all EQUIPMENT and SYSTEMS

Guidance and manufacturer's declaration-electromagnetic immunity			
<i>Vital signs monitor</i> is intended for use in the electromagnetic environment specified below. The user of <i>VITAL SIGNS MONITOR</i> should make sure that it is used in such environment.			
Immunity test	EN 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge(ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	2kV for power supply lines 1 kV for input/output lines	2kV for power supply lines 1kV for input/output lines	If display data reading is not stable, should use spare stable-pressure supply and try to make the functional ground line and supply cable collateral line with close quarters.
Surge IEC 61000-4-5	1 kV line(s) to line(s) 2 kV line(s) to earth	1 kV line(s) to line(s) 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % <i>UT</i> (>95 % dip in <i>UT</i>) for 0,5 cycle 40 % <i>UT</i> (60 % dip in <i>UT</i>) for 5 cycles 70 % <i>UT</i> (30 % dip in <i>UT</i>) for 25 cycles <5 % <i>UT</i> (>95 % dip in <i>UT</i>) for 5 sec	<5 % <i>UT</i> (>95 % dip in <i>UT</i>) for 0,5 cycle 40 % <i>UT</i> (60 % dip in <i>UT</i>) for 5 cycles 70 % <i>UT</i> (30 % dip in <i>UT</i>) for 25 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the <i>VITAL SIGNS MONITOR</i> requires continued operation during power mains interruptions, it is recommended that the <i>VITAL SIGNS MONITOR</i> be powered from an uninterruptible power supply or a battery.

Power frequency (50Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
---	------	------	---

**Guidance and manufacturer’s declaration-electromagnetic immunity
For EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING**

Guidance and manufacturer’s declaration-electromagnetic immunity			
<i>Vital signs monitor</i> is intended for use in the electromagnetic environment specified below. The user of <i>Vital signs monitor</i> should make sure that it is used in such environment.			
Immunity test	EN 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the <i>Vital signs monitor</i> , including cables. The recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \left[\frac{3.5}{V1} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E1} \right] \sqrt{P} \quad \mathbf{80\ MHz\ to\ 800\ MHz}$ $d = \left[\frac{7}{E1} \right] \sqrt{P} \quad \mathbf{800\ MHz\ to\ 2.5\ GHz}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

			<p>Where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitter as determined by an electromagnetic site survey ^a should be less than the compliance level in each frequency range ^b. Interference may occur in the vicinity of equipment marked with the following symbol:</p>
<p>ATTENTION 1 At 80 MHz. the higher frequency range applies.</p>			
<p>ATTENTION 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p>a. Field strengths from fixed transmitters, such as base stations for radio(cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters. an electromagnetic site survey should be considered. If the measured field strength in the location in which the VITAL SIGNS MONITOR is used exceeds the applicable RF compliance level above, the VITAL SIGNS MONITOR should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the VITAL SIGNS MONITOR.</p> <p>b. Over the frequency range 150 kHz to 80MHz. Field strengths should be less than 1 V/m.</p>			

Recommended separation distances between portable and mobile RF RF communications equipment and the EQUIPMENT or SYSTEM For EQUIPMENT or SYSTEM that are not LIFE-SUPPORTING

<p>Recommended separation distances between Portable and mobile RF communications equipment and the VITAL SIGNS MONITOR</p>	
<p><i>Vital signs monitor</i> is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user of <i>Vital signs monitor</i> can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the <i>Vital signs monitor</i> as recommended below according to the maximum output power of the communications equipment.</p>	
<p>Rated</p>	<p>Separation distance according to frequency of transmitter(m)</p>

maximum output power of transmitter (W)	150 kHz to 80 MHz $d = \left[\frac{3.5}{V1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{E1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E1} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.75
1	1.17	1.17	2.33
10	3.69	3.69	7.37
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter. Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

ATTENTION1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

ATTENTION2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. However, because of the proliferation of radio-frequency transmitting equipment and other sources of electrical noise in the health-care and home environments (for example, cellular phones, mobile two-way radios, electrical appliances), it is possible that high levels of such interference due to close proximity or strength of a source, may result in disruption of performance of this device.

Vital signs monitor generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with these instructions, may cause harmful interference with other devices in the vicinity. Disruption may be evidenced by erratic readings, cessation of operation, or other incorrect functioning. If this occurs, the site of use should be surveyed to determine the source of this disruption. Take below actions to eliminate the source:

- Turn the equipment in the vicinity off and on to isolate the offending equipment.
- Reorient or relocate the other receiving devices.
- Increase the separation between the interfering equipment and VITAL SIGNS MONITOR.

If assistance is required, contact Technical Services Department or your local representatives.

Chapter 9: Specification

9.1 Equipment Classification

Type of Protection Against Electric shock:	Internally Powered
Mode of operation:	Continuous
Degree of Protection Against Ingress of Liquids:	IPX1, drip proof
Degree of Mobility:	Portable
Degree of Protection Against Electric Shock:	Type BF
Safety Requirements:	EN60601-1:1990

9.2 Parameter

Parameter	Range	Accuracy	Resolution
SpO₂	0-100%	±2 at 70 - 100% <70%, undefined	1%
Pulse Rate	25-250bpm	±2 at 25 – 250bpm	1bpm
EtCO₂	0-20% (0-150mmHg)	±2 mm Hg @ < 5.0% CO ₂ (at BTPS) < 10% of reading @ > 5.0% CO ₂ (at BTPS)	0.1 mmHg
RR	0-150 (BPM)	±1bpm	1bpm

9.6 Default Settings of Alarms Limits

	High Alarm Limits			Low Alarm Limits		
	Adult	Pediatric	Neonatal	Adult	Pediatric	Neonatal
SpO₂	99	99	99	92	92	92
PR	100	110	120	60	70	80
ETCO₂	40	40	40	30	30	30
RR	30	30	30	8	8	8

9.7 Power Requirements

AC power supply 100-240VAC, 47-63Hz

Or four standard “AA” alkaline or Ni-MH cells (IEC Type LR6)

9.8 Battery Life

Alkaline Cells: approximately 6 hours

9.9 Dimensions

Width	75mm (2.95 inches)
Height	135mm (5.31 inches)
Depth	28mm (1.10 inches)
Weight	258grams (9.10 ounces) with batteries

9.10 Environmental Specification

Operating Temp.	0 to 45°C
Storage Temp.	-40 to +55°C
Relative Humidity	30 to 95% (operating) 10 to 95% (storage) 0 to 95% (shipping)